



الكلية: الطب

القسم او الفرع: الاحياء المجهرية

المرحلة: الثالثة

أستاذ المادة: د. ميساء إبراهيم محمد

اسم المادة باللغة العربية: احياء مجهرية

اسم المادة باللغة الإنكليزية: **Microbiology**

اسم المحاضرة الرابعة باللغة العربية: الاختبارات المناعية

اسم المحاضرة الرابعة باللغة الإنكليزية: **Practical Lecture on Immunology Techniques**

**PRACTICAL LECTURE ON
IMMUNOLOGY TECHNIQUES**

Dr- Maysaa Ibrahim

Procedure:

1. PRECIPITATION REACTIONS

Precipitation occurs when a soluble antigen reacts with its specific antibody to form an insoluble complex visible as a precipitate

A. Ring Test:

B. Principle:

C. When antigen and antibody solutions are layered, a visible precipitate forms at their interface.

D. Materials:

E. Antigen solution

F. Antiserum (antibody)

G. Narrow test tubes or small tubes

H. Micropipettes

1-Place 1 mL of antiserum in a test tube.

2- Carefully overlay 1 mL of antigen solution without mixing.

3- Allow the tube to stand at room temperature undisturbed.

Observation:

A thin white ring of precipitate forms at the junction

if the antigen and antibody are specific.

A. Slide Agglutination

1- Precipitation Reactions

B. Single Radial Immunodiffusion (Mancini Method)

C- Double Immunodiffusion (Ouchterlony Test)

D. Immunoelectrophoresis

2. Agglutination Reactions

Agglutination occurs when insoluble antigens (cells, particles) react with specific antibodies, causing clumping.

Principle:

Particulate antigens (e.g., bacteria) clump when mixed with specific antibodies on a slide.

Materials:

- Glass slides**
- Bacterial suspension**
- Antiserum**
- Mixing sticks**

Procedure:

- 1.Place a drop of antigen suspension on a clean slide.**
- 2.Place a drop of antiserum next to it.**
- 3.Mix gently with a stick.**
- 4.Rock the slide and observe.**

Observation:

- Visible clumping within minutes indicates a positive result.

B. Tube Agglutination (Widal Test Example)

Principle:

Serial dilutions of serum are mixed with antigen suspensions in tubes to detect and measure antibody titers.

Materials:

- Test tubes
- Patient serum
- Known antigen suspensions
- Normal saline

Procedure:

- 1.Make serial two-fold dilutions of serum in tubes.
- 2.Add an equal volume of antigen suspension to each tube.
- 3.Incubate at 37°C for 2–4 hours or overnight.

Observation:

- Clumping/sedimentation in tubes indicates agglutination.
- Titer = highest dilution showing visible agglutination.

3- Complement Fixation Test (CFT)

Principle:

Complement is fixed if antigen-antibody reaction occurs. If complement remains free, it causes hemolysis of indicator red cells.

Materials:

- Complement source (e.g., guinea pig serum)

- Antigen
- Patient serum (antibody)
- Sensitized sheep RBCs (indicator system)

Procedure:

- 1.Incubate antigen with patient serum and complement.
- 2.Add sensitized sheep RBCs.
- 3.Observe hemolysis.

Interpretation:

- No hemolysis = Positive test (complement fixed, antigen-antibody present).
- Hemolysis = Negative test (no antigen-antibody reaction).

Thank you