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**Serial dilution is the stepwise dilution of a substance in solution. Usually, the dilution factor at each step is constant.**

**A ten-fold serial dilution could be 1 M, 0.1 M, 0.01 M, 0.001 M ... Serial dilutions are used to accurately create highly diluted solutions as well as solutions for experiments resulting in concentration curves with a logarithmic scale.**

**Building of standard curve:**

**After taking the absorbance by UV spectroscopy device, make (x , y) points on graph paper representing (Conc. , ABS).**

#	Conc.	ABS
1		
2		
3		
4		
5		

**Draw a straight line that passes among the points without deviation (not necessarily to hit all the points).**

**The standard curve is used to determine the concentration of any sample containing the same compound used in serial dilution, by obtaining the absorbance of that sample.**

**This can be done either by determination on plotting, or using the straight-line equation.**

**Least square method:**

$$\text{slope} = \frac{\sum[(x - x') * (y - y')]}{\sum[(x - x')^2]}$$

**x' is the average of x values**

**y' is the average of y values**

**UV spectrophotometer:**

**Cuvette**

**Blank**

**Lambda max ( $\lambda_{\text{max}}$ )**

**Absorbance**

**Example:**

**The following results are documented after serial dilution of salicylic acid stock solution (200mg/100ml).**

**Build a standard curve, then find the conc. of a sample that has absorbance of 0.54 Au.**

#	Conc. ( $\mu\text{g/ml}$ )	ABS (Au)
1	200	0.75
2	20	0.35
3	2	0.2
4	0.2	0.09
5	0.02	0.01