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Chlorobutanol synthesis

1-Introduction

Chlorobutanol is 1,1,1-trichloro-2-methyl-2-propanol (it is a tertiary alcohol)

It has been used in therapeutic for a variety of purposes:

1-Bacterostatic (used as preservative) in many injectable, ophthalmic & intranasal preparations.

2-Sedative-hypnotic in the past.


3-Local anesthetic, in painful I.M injectable, topical preparations & as a dental analgesic

2-Physical properties:

1-White crystalline substance exists in two forms (anhydrous & hydrated)

2-Characteristic camphor-like odour & taste


3-Slightly soluble in H₂O (1/125) & freely soluble in alcohol (1/1). It is more soluble in boiling H₂O but hydrolysis at such temp. will take place readily.

Therefore recrystallized from H₂O-alcohol mixture.


3-Preparation :

It is prepared from acetone & chloroform in alcoholic KOH (to speed the reaction)

Synthesis Chlorobutanul from Furan

4-Chloro-1-butanol is obtained simply by leading [hydrogen chloride](#) into boiling tetrahydrofuran until the temperature of the boiling mixture reaches 103.5-105.5°C (ca. 5 hours). 4-Chloro-1-butanol is isolated in 54-57% yield by distillation in a vacuum, at 15 mm, much hydrochloric acid being at first evolve. 4-Chloro-1-butanol 84-85 °C/16 mmHg.

4-Procedure:

1-Mix 33 ml acetone (25 gm) in dry conical flask with 7 ml chloroform (10 gm)

2-Cool the mixture

3-Prepare alcoholic KOH (1.75 gm KOH / 12.5 ml alcohol) & put in separatory funnel

4-Add this solution (alc. KOH) drop wise from the separatory funnel on the previous solution within a period of 15 min.

5-Filter the ppt. KCl

6-Evaporate the filtrate on H₂O bath